

# ABDOMINAL BINDERS

## Deluxe Abdominal Binders 1231, 1241, 1251

Designed for support after surgery and in specific types of hernias. May also be used for treatment of back strain. Acts as secondary dressing. Velcro® closure simplifies application for professional and patient. Double darts facilitate contour fit. Sized from small through extra-extra-large. See size measurements below.

- 1231-0**                    9" Width Abdominal Binders
- 1241-0**                    10" Width Abdominal Binders
- 1251-0**                    12" Width Abdominal Binders



## Abdominal Binder 1241-1

Manufactured of knitted elastic and closed with Velcro. This economical binder has double darts in the back. 10" wide elastic is sized by hip measurements from small through extra-extra-large.

- 1241-1**                    10" Elastic Binder

*For all of the sized binders listed above, measure hip area and order as follows:  
Small 24" - 30" - Medium 30" - 36" - Large 36" - 42"  
Extra-Large 42" - 48" - Extra-Extra-Large 48" - 54"*

## Deluxe Segmented Abdominal Binders 1423

3 or 4 each 3" plushbacked elastic panels are stitched together for form fitting excellence. Strong, soft elastic provides the necessary support and Velcro closure provides for size adjustment. #1423-00, the 9" wide binder, has 3 elastic panels as shown. #1423-20, the 12" wide binder has 4 panels. One regular universal size fits both male and female hip measurements of 30" through 45". An extra-large universal is also available for 46" through 62" hips.

- 1423-00**                    9" Deluxe Binder, Universal
- 1423-04**                    9" Deluxe Binder, Extra-large Universal
- 1423-20**                    12" Deluxe Binder, Universal
- 1423-24**                    12" Deluxe Binder, Extra-large Universal



## Segmented Abdominal Binders 2600

The 3" wide, knitted elastic segmented sections provide form fitting comfort. Velcro closure. The 9" binder has 3 panels and the 12" binder has 4 panels. Each width is available in a regular-universal size which fits males or females with hip measurements of 30" through 45", an extra-large universal size for hip measurements of 46" through 62" and a two "X" size for 63" through 75".

- 2600-00**                    9" Elastic Binder, Universal
- 2600-04**                    9" Elastic Binder, Extra-large Universal
- 2600-05**                    9" Elastic Binder, Extra-extra-large Universal
- 2600-20**                    12" Elastic Binder, Universal
- 2600-24**                    12" Elastic Binder, Extra-large Universal
- 2600-25**                    12" Elastic Binder, Extra-extra-large Universal